**CRYOcheck™** **IVD**  
**WEAK LUPUS POSITIVE CONTROL**  
PLASMA DE CONTROLE  
LUPUS POSITIF FAIBLE

**CRYOcheck™**  
**Hemostasis Control Plasmas**  
**WEAK LUPUS POSITIVE CONTROL**

**Intended Use**

CRYOcheck Weak Lupus Positive Control is prepared from human source plasma and is recommended for use as a positive control in assays for lupus anticoagulant.

**Summary and Principle**

An association between circulating anticoagulants and systemic lupus erythematosus (SLE) was first described in 1952<sup>1</sup>. In 1972 the term “lupus anticoagulant” (LA) was used to describe these non-specific, circulating inhibitors<sup>2</sup>. The designation LA is still used today although it is now evident the majority of these patients do not suffer from SLE<sup>3</sup>. The presence of LA has been increasingly associated with a variety of disorders such as unexplained thrombosis and recurrent fetal loss<sup>4,5</sup>.

LA are immunoglobulins, usually of the IgG, IgM or IgA classes, which are directed against anionic phospholipids or phospholipid-plasma protein complexes causing interference with *in vitro* phospholipid dependent coagulation tests. LA are clinically distinct members of a broader group of antiphospholipid antibodies (APA) characterized by antigenic protein targets such as β<sub>2</sub> glycoprotein 1 and prothrombin. APA include LA as well as anticardiolipin, antiphosphatidylserine, and antiphosphatidylethanolamine antibodies<sup>6,7,8</sup>.

The SSC Subcommittee has specified diagnostic criteria for the detection of LA<sup>9</sup>.

**Reagents**

CRYOcheck Weak Lupus Positive Control contains citrated human plasma collected from donors that have tested positive in accordance with the revised criteria of the SSC Subcommittee for the Standardization of Lupus Anticoagulants<sup>9</sup>. Source plasmas are processed in a manner that yields platelet-poor plasmas. Plasma is then buffered, aliquoted and frozen at very low temperatures.

 *All blood products should be treated as potentially infectious. Source material from which this product was derived was found to be negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. Accordingly, these human blood based products should be handled and discarded as recommended for any potentially infectious human specimen<sup>10</sup>.*

**Storage and Handling**

When stored at -40 to -80 °C, CRYOcheck Weak Lupus Positive Control is stable to the end of the month indicated on the product packaging.

Thaw each vial at 37 °C (± 1 °C) in a waterbath. **The use of a dry bath or heating block for thawing is not recommended.** Thawing times are important and should be strictly adhered to. The use of a timer is recommended. Refer to the Thawing Table for recommended thawing times based on aliquot size. Allow thawed plasma to acclimate to room temperature (18 to 25 °C) and invert gently prior to use.

Thawing Table	
Aliquot Size	37 °C (± 1 °C) Waterbath
0.5 mL	3 minutes
1.0 mL	4 minutes

CRYOcheck Weak Lupus Positive Control may be used for up to eight hours after thawing, if capped in the original vial and maintained at 2 to 8 °C. Allow refrigerated plasma to acclimate to room temperature (18 to 25 °C) and invert gently prior to use. **Thawed material should be discarded after eight hours and should not be refrozen.**

Availability

Product	Catalog #	Format
Weak Lupus Positive Control	CCWLP-05	25 vials x 0.5 mL
	CCWLP-10	25 vials x 1.0 mL

**Instruments**

Each lab should prepare the local instrument in accordance with the manufacturer’s instructions for use.

**Procedure**

After thawing and preparing CRYOcheck Weak Lupus Positive Control, use in accordance with established laboratory procedures for the quality control of assays for LA.

Materials Provided

- CRYOcheck Weak Lupus Positive Control

Materials Required but not Provided

- Waterbath capable of maintaining 37 °C (± 1 °C)
- Assay reagents
- Coagulation instrument or assay system
- Sample cups
- Plastic disposable pipettes
- Volumetric pipette

Quality Control

Each laboratory should establish its own quality control (QC) ranges using acceptable statistical methods. These QC ranges may then be used to monitor and validate the integrity of the test system<sup>11</sup>. For all coagulation tests, the laboratory must include at least two levels of control for every eight hours of operation and any time a change in reagents occurs<sup>12</sup>.

**Results**

Control results should fall within the laboratory’s established QC ranges provided the integrity of the test system has not been compromised.

**Limitations of the Procedure**

When proper control values are not obtained, assessment of each component of the test system including reagents, control plasmas, instrumentation and operator technique must be undertaken in order to ascertain that all other components are functioning properly.

**Expected Values**

Quality control testing yielded the following results with a representative lot of CRYOcheck Weak Lupus Positive Control. Results may vary depending on instrument, method and technique used.

Lupus Anticoagulant Assays	Result	Interpretation
APTT (Lupus Sensitive)	69.3 sec	Positive
APTT (1:1 Normal Mix)	46.2 sec	Positive
APTT (1:1 Saline Mix)	67.7 sec	Positive
Kaolin Clotting Time (KCT)	166 KCT units	Positive
dilute Russell’s Viper Venom Time (dRVVT)	1.3 ratio	Positive
Hexagonal Lupus Anticoagulant	14.9 sec	Positive
Platelet Neutralization Procedure (PNP)	16.0 sec	Positive

**Performance Characteristics**

When used in accordance with established testing procedures for LA<sup>9</sup>, CRYOcheck Weak Lupus Positive Control demonstrates the following characteristics:

- CRYOcheck Weak Lupus Positive Control shows prolongation in phospholipid-dependent clotting tests;
- The clotting time of a mixture of CRYOcheck Weak Lupus Positive Control and pooled normal plasma is significantly longer (≥3 SD) than that of pooled normal plasma and non-LA patient plasma
- CRYOcheck Weak Lupus Positive Control is corrected by the addition of lysed, washed platelets or phospholipid liposomes containing phosphatidylserine or hexagonal phase phospholipids;
- CRYOcheck Weak Lupus Positive Control is non-specific for any individual clotting factor.

Each laboratory should establish their own diagnostic criteria for LA assays. CRYOcheck Weak Lupus Positive Control is subject to the limitations of the assay system in use.

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**Bibliography / Bibliographie**

- Conley CL, Hartmann RC. A hemorrhagic disorder caused by circulating anticoagulant in patients with disseminated lupus erythematosus. J Clin Invest 1952; 31:621-622.
- Feinstein DI, Rapaport SI. Acquired inhibitors of blood coagulation. Progress in hemostasis and thrombosis 1972; 1:75-95.
- Schleider MA, Nachman RL, Jaffe EA, Coleman M. A clinical study of the lupus anticoagulant. Blood 1976; 48:499-509.
- Bowie EJW, Thompson JH, Jr., Pascuzzi CA, Owen CA, Jr. Thrombosis in systemic lupus erythematosus despite circulating anticoagulants. J Lab Clin Med 1963; 62:416-430.
- Nilsson IM, Astedt B, Hedner U, Berezin D. Intrauterine death and circulating anticoagulant, “Antithromboplastin”. Acta Med Scand 1975; 197(3):153-159.
- Thiagarajan P, Shapiro SS, DeMarco L. Monoclonal immunoglobulin M lambda coagulation inhibitor with phospholipid specificity: Mechanism of a lupus anticoagulant. J Clin Invest 1980; 66(3):397-405.
- Triplet DA. Coagulation assays for the lupus anticoagulant: Review and critique of current methodology. Stroke 1992; 23(I); I-11 – I-14.
- Permpikul P, Rao LVM, Rapaport SI. Functional and binding studies of the roles of prothrombin and β<sub>2</sub> glycoprotein I in the expression of lupus anticoagulant activity. Blood 1994; 83(10):2878-2892.
- Pengo V, Tripodi A, Reber G, Rand JH, Ortel TL, Galli M, De Groot PG; Update of the guidelines for lupus anticoagulant detection. Subcommittee on Lupus Anticoagulant/Antiphospholipid Antibody of the Scientific and Standardisation Committee of the International Society on Thrombosis and Haemostasis. J Thromb Haemost. 2009; 7(10):1737-1740.
- Biosafety in Microbiological and Biomedical Laboratories 5<sup>th</sup> ed. Centers for Disease Control and Prevention / National Institutes of Health, 2009.
- Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989. p. 166-171.
- CLIA 2004 – Code of Federal Regulations, 42CFR493.1269, 2004.

 European Authorized Representative (Regulatory affairs only)  
Emergo Europe—Prinsessegracht 20, 2514 AP The Hague, The Netherlands

*Symbols used / Symboles utilisés*

09.60.00016

Rev. 08 May / mai 2017

*In vitro* diagnostic medical device

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