## CRYOcheck™ IVD NORMAL REFERENCE PLASMA

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## CRYOcheck™ Calibration Plasma NORMAL REFERENCE PLASMA

*A Member of our Gold Standard Family of Assayed Plasmas*

#### Intended Use

CRYO*check* Normal Reference Plasma is recommended for use in the *in vitro* quantification of hemostatic parameters in human plasma.

#### Summary and Principle

The use of reference plasma is widely recommended for the quantification of hemostatic parameters in human plasma<sup>1-4</sup>. These procedures are performed as part of the laboratory evaluation of patients with coagulation disorders and require the construction of reference or dose response curves from which quantitative measures of individual analytes can then be determined. The World Health Organization (WHO) has established a series of international standards for this purpose in an attempt to standardize these procedures.

#### Reagents

CRYO*check* Normal Reference Plasma consists of normal citrated human plasma collected from a minimum of 20 carefully screened donors. The plasma pool is buffered using HEPES buffer, aliquoted and rapidly frozen. CRYO*check* Normal Reference Plasma exhibits hemostatic parameters which are representative of a normal population. Each lot number is assayed using international reference standards (where available) and values for hemostatic parameters are assigned. *Refer to the ASSAY CERTIFICATE for the assigned values specific to each lot number.*

 *All blood products should be treated as potentially infectious. Source material from which this product was derived was found to be negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. Accordingly, these human blood based products should be handled and discarded as recommended for any potentially infectious human specimen*<sup>5</sup>.

#### Storage and Handling

When stored at -40 to -80 °C, CRYO*check* Normal Reference Plasma is stable to the end of the month indicated on the product packaging.

Thaw each vial at 37 °C (± 1 °C) in a waterbath. **The use of a dry bath or heating block for thawing is not recommended.** Thawing times are important and should be strictly adhered to. The use of a timer is recommended. Refer to the Thawing Table for recommended thawing times based on aliquot size. Allow thawed plasma to acclimate to room temperature (18 to 25 °C) and invert gently prior to use.

Thawing Table	
Aliquot Size	37 °C (± 1 °C) Waterbath
0.5 mL	3 minutes
1.0 mL	4 minutes

CRYO*check* Normal Reference Plasma may be used for up to eight hours after thawing, if capped in the original vial and maintained at 2 to 8 °C. Allow refrigerated plasma to acclimate to room temperature (18 to 25 °C) and invert gently prior to use. **Thawed material should be discarded after eight hours and should not be refrozen.**

#### Availability

Product	Catalog #	Format
Normal Reference Plasma	CCNRP-05	25 vials x 0.5 mL
	CCNRP-10	25 vials x 1.0 mL

#### Instruments

Each lab should prepare the local instrument in accordance with the manufacturer’s instructions for use.

#### Procedure

After thawing and preparing CRYO*check* Normal Reference Plasma, use in accordance with established laboratory procedures for the quantitative assessment of hemostatic parameters.

#### Materials Provided

- CRYO*check* Normal Reference Plasma

#### Materials Required but not Provided

- Waterbath capable of maintaining 37 °C (± 1 °C)
- Timer
- Assay reagents
- Coagulation instrument or assay system
- Quality control material (e.g. CRYO*check* Reference Control Normal, CRYO*check* Abnormal 1 Reference Control, CRYO*check* Abnormal 2 Reference Control)
- Sample cups
- Volumetric pipette
- Plastic disposable pipettes

#### Quality Control

Each laboratory should establish its own quality control (QC) ranges using acceptable statistical methods. These QC ranges may then be used to monitor and validate the integrity of the test system<sup>6</sup>. For all coagulation tests, the laboratory must include at least two levels of control for every eight hours of operation and any time a change in reagents occurs<sup>7</sup>.

#### Results

Control results should fall within the laboratory’s established QC ranges provided the integrity of the test system has not been compromised.

#### Limitations of the Procedure

When proper control values are not obtained, assessment of each component of the test system including reagents, control plasmas, instrumentation and operator technique must be undertaken in order to ascertain that all other components are functioning properly.

#### Expected Values

Refer to the **ASSAY CERTIFICATE** for the expected values specific to each lot number of CRYO*check* Normal Reference Plasma.

#### Performance Characteristics

Values that have been assigned to CRYO*check* Normal Reference Plasma have been determined in accordance with accepted clinical laboratory procedures using the methods indicated on the **ASSAY CERTIFICATE**. All components in each individual system should be assessed to determine their effect on the reproducibility and accuracy of expected values. When used properly, CRYO*check* Normal Reference Plasma is subject to the limitations of the assay system in use.

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## PrecisionBioLogic

**Bibliography / Bibliographie**

1. Dombrose FA, Barnes CC Jr, Gaynor JJ, Elston RC. A lyophilized human reference plasma for coagulation factors. Evidence for stability of factors I, II, V, and VII through XII. Am J Clin Pathol. 1982; 77(1):32-45.
2. Burgi W, Schnell, E. Artificial control materials: Coagulation. In: Rosalki, SB, editor. New approaches to laboratory medicine. Darmstadt: G-I-T Verlag Ernst Giebeler; 1981. p. 57-65.
3. Thelin GM. Preparation and standardization of a stable AHF plasma. Thromb Diath Haemorrh 1968; 19(3):423-429.
4. Dombrose FA, Barnes CC. A standard reference plasma for coagulation assays. In: Triplett DA, editor: Standardization of coagulation assays: An overview. Illinois: ASCP Press; 1982. p. 223-234.
5. Biosafety in Microbiological and Biomedical Laboratories 5<sup>th</sup> ed. Centers for Disease Control and Prevention / National Institutes of Health, 2009.
6. Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989. p. 166-171.
7. CLIA 2004 – Code of Federal Regulations, 42CFR493.1269, 2004.

 European Authorized Representative (Regulatory affairs only)  
Emergo Europe—Prinsessegracht 20, 2514 AP The Hague, The Netherlands

*Symbols used / Symboles utilisés*

09.60.00013

Rev. 11 May / mai 2017

*In vitro* diagnostic medical device

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