Quantitative Protein S Clotting Assay
CLOT S™

Intended Use

CRYOcheck Clot S is a clot-based assay intended for the quantitative determination of protein S activity in citrated human plasma.

Summary and Principle

Protein S is a vitamin K-dependent glycoprotein with a molecular weight of ~70,000 Da. The mature protein circulates in plasma at a total concentration of 25 µg/mL, 40% free and 60% reversibly bound to C4b-binding protein¹. The free form of protein S acts as a cofactor to enhance the anticoagulant property of activated protein C, a potent inhibitor of procoagulant factors Va and VIIIa².

Protein S deficiency has both congenital and acquired etiologies of clinical interest. Congenital deficiencies are commonly associated with an increased risk of venous thromboembolic disease^{3,4,5}, while acquired abnormalities are associated with anti-vitamin K therapy⁶, liver disease⁶, disseminated intravascular coagulation (DIC)⁶, oral contraceptives⁷, oestrogen therapy⁸, acute phase inflammatory responses⁹, pregnancy¹⁰, and newborns¹¹.

CRYOcheck Clot S initiates the common pathway of coagulation in plasma using a Russell's viper venom (RVV-X) reagent to convert factor X to Xa in the presence of activated protein C (APC), bypassing all factors above the common pathway¹². When mixed with protein S deficient plasma, samples from patients with a protein S deficiency or dysfunction will have shortened CRYOcheck Clot S clotting times relative to samples with normal levels of functional protein S. The clotting time is proportional to the amount of functional protein S in the patient's plasma and this can be quantified using a calibration curve.

Reagents

Protein S Deficient Plasma (PS Deficient)

Contains citrated pooled human plasma (depleted of protein S by immunoabsorption), buffers and stabilizers.

Clot S Activator (Activator)

Contains activated protein C, Russell's viper venom, heparin neutralizing agents, buffers and stabilizers.

C & S Diluent

Available separately from Precision BioLogic (catalog # CSD).

All blood products should be treated as potentially infectious. Source material from which this product was derived was found to be negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. Accordingly, these human blood based products should be handled and discarded as recommended for any potentially infectious human specimen¹³.

Storage, Preparation and Handling

When stored at -70 °C or below, CRYOcheck Clot S is stable to the end of the month indicated on the product packaging.

Thaw one vial each of **PS Deficient** and **Activator** at 37 °C (± 1 °C) in a waterbath using the "floatie" thawing device (provided separately). Thawing times are important and should be strictly adhered to. The use of a timer is recommended. Refer to the Thawing Table for recommended thawing times according to format.

Thawing Table	
Aliquot Size	37 °C (± 1 °C) Waterbath
3.0 mL	6 minutes
1.5 mL	5 minutes

Immediately after thawing:

- vortex or vigorously mix **PS Deficient** only for 5 - 10 seconds
- gently swirl **Activator** for 5 - 10 seconds
- allow uncapped reagents to acclimate on-board the instrument for **30 minutes**

- Gently mix **PS Deficient** prior to performing tests and every two hours thereafter to prevent reagent sedimentation. For detailed mixing instructions, please consult the protocol for your automated coagulation analyzer.

For manual methods, swirl **PS Deficient** prior to performing tests and every two hours thereafter.

CRYOcheck Clot S may be used for up to five hours after preparation. When not in use, CRYOcheck Clot S reagents should be capped in the original vials and maintained at 2 to 8 °C. Before use, repeat steps 1 to 4.

NB: CRYOcheck Clot S Protein S Deficient Plasma and Clot S Activator are lot specific and should not be interchanged with other lot numbers.

Availability

Product	Catalog #	Format	# of Tests
Clot S	CCS-30	PS Deficient 5 x 3.0 mL Activator 5 x 3.0 mL	300
	CCS-15	PS Deficient 5 x 1.5 mL Activator 5 x 1.5 mL	150

Instruments

Each lab should prepare the local instrument in accordance with the manufacturer's instructions for use.

Procedure

Materials Provided

- Protein S deficient Plasma (**PS deficient**)
- Clot S activator (**Activator**)

Materials Required but not Provided

- C & S Diluent
- 0.025 M CaCl₂
- Waterbath capable of maintaining 37 °C (± 1 °C)
- Floatie for thawing vials in waterbath
- Coagulation instrument or assay system
- Calibration plasma (e.g. CRYOcheck Normal Reference Plasma)
- Quality control material (e.g. CRYOcheck Reference Control Normal, CRYOcheck Abnormal 1 Reference Control, CRYOcheck Abnormal 2 Reference Control)
- Graph paper
- Plastic test tubes (e.g. 12 x 75 mm)
- Coagulation reaction cuvettes
- Plastic disposable pipettes
- Volumetric pipette
- Timer

Specimen Collection and Preparation

Patient samples should be collected into 105 - 109 mmol/L sodium citrate dihydrate anticoagulant (3.2%) in a ratio of 9 parts blood to 1 part anticoagulant. Patient plasma is derived by centrifugation at 1,500 x g for 15 minutes in order to achieve platelet-poor plasma (<10,000 platelets/µL) and should be tested within four hours of collection when maintained at 2 to 4 °C. If samples are not to be tested within four hours then plasma should be removed from the cells and frozen at -20 °C for up to two weeks or -70 °C for up to six months in accordance with the Clinical Laboratory Standards Institute (CLSI) guidelines¹⁴.

Assay Procedure

- Prepare CRYOcheck Clot S reagents according to Storage, Preparation and Handling instructions above.
- Prepare instrument according to the manufacturer's instructions for use.
- Prepare a 1:10 dilution of test plasma (i.e. patient, calibrator or control) in C & S Diluent (**do not substitute distilled water or other buffers for C & S Diluent**).
- To a coagulation reaction cuvette, immediately add 50 µL of the diluted test plasma, 50 µL of PS Deficient and 50 µL of Activator.
- Mix and incubate at 37 °C for three minutes.
- Add 50 µL 0.025 M CaCl₂ and immediately initiate timer.
- Record clotting time in seconds.

Assay Calibration


- Prepare CRYOcheck Clot S reagents according to Storage, Preparation and Handling instructions above.
- Prepare calibration plasma according to manufacturer's directions.
- Prepare serial dilutions of calibration plasma from 1:7 to 1:80 in C & S Diluent according to the following table:

Tube No.	C & S Diluent (mL)	Calibration Plasma (mL)	Dilution	% Protein S
1	0.6	0.1	1:7	143
2	1.8	0.2	1:10	100
3	1.0	1.0 of Tube No. 2	1:20	50
4	1.0	1.0 of Tube No. 3	1:40	25
5	1.0	1.0 of Tube No. 4	1:80	12.5
6	1.0	0	n/a	0

Note: This is an example only of a serial dilution profile prepared using calibration plasma with a protein S level of 100%. Always be sure to utilize the lot-specific % protein S level of the calibration plasma in use. If using CRYOcheck Normal Reference Plasma, refer to the lot-specific Assay Certificate.

- To a coagulation reaction cuvette, add 50 µL from Tube 1, 50 µL of PS Deficient, and 50 µL of Activator. Mix and incubate at 37 °C for three minutes.
- Add 50 µL 0.025 M CaCl₂ and immediately initiate timer. Record clotting time in seconds.
- Repeat steps 4 and 5 for Tubes 2 through 6.
- Plot clotting times in seconds (y-axis) vs. % of protein S activity (x-axis).
- Construct a standard curve and derive % protein S values (see Example Only: CRYOcheck Clot S Calibration Curve). It is recommended that samples exceeding 140% protein S be diluted at 1:20 and re-tested. Samples below 10% protein S should be diluted at 1:5 and re-tested.

Example Only: CRYOcheck Clot S Calibration Curve

Quality Control

Each laboratory should establish its own quality control (QC) ranges using acceptable statistical methods. These QC ranges may then be used to monitor and validate the integrity of the testing system¹⁵. For all coagulation tests, the laboratory must include at least two levels of control for every eight hours of operation and any time a change in reagents occurs¹⁶.

Results

Results are expressed as a percentage of normal protein S activity by comparison with a known standard or calibration plasma. Protein S values recovered below the laboratory established normal range may be indicative of a protein S deficiency (congenital or acquired). Each laboratory should establish its own normal reference range for protein S activity in accordance with CLSI guidelines¹⁷.

Limitations of the Procedure

Factor VIII:c Interference: CRYOcheck Clot S is unaffected by factor VIII:c activity up to 600%.

Heparin Interference: CRYOcheck Clot S is unaffected by unfractionated heparin (UFH) or by low molecular weight heparin (LMWH) up to 1.0 IU/mL.

Direct Thrombin Inhibitors: CRYOcheck Clot S may be affected by hirudin and other direct thrombin inhibitors, resulting in falsely elevated protein S activity levels.

Lupus Anticoagulant: Lupus anticoagulants (LA) present in test plasma may affect CRYOcheck Clot S results.

Activated Protein C Resistance: CRYOcheck Clot S results may be affected by samples from patients with the factor V_{Leiden} mutation.

Expected Values

A normal population study was performed on 100 healthy adults. A mean protein S level of 90.3% with a two standard deviation (SD) range of 46.7% – 133.9% was recovered. It is recommended that each laboratory establish its own normal population range.

Performance Characteristics

Reportable Range: A reportable range for CRYOcheck Clot S of 10 to 140% protein S activity was obtained on a STA-R® instrument (Diagnostica Stago). Reportable range may vary depending on instrumentation used.

Precision: Intra-assay reproducibility was assessed by testing one normal and one abnormal plasma (with reduced % protein S) 20 times each. To evaluate inter-assay precision, one normal sample and one abnormal sample were tested over multiple days, using multiple calibration curves. Mean, SD, and percent coefficient of variation (%CV) were as follows:

Test Sample	Intra-Assay Precision			Inter-Assay Precision		
	Mean	SD	%CV	Mean	SD	%CV
Normal	90.5	2.1	2.3	101.6	6.8	6.7
Abnormal	26.3	1.9	7.2	34.9	2.8	7.9

Correlation: CRYOcheck Clot S was compared to another commercially available clot-based protein S test, using clinical samples from the assay target population.


Figure 1: Correlation of protein S values determined on 100 samples from the target population.

Bibliography / Bibliographie

- DiScipio RG, Davie EW. Characterization of protein S, a gamma-carboxyglutamic acid containing protein from bovine and human plasma. Biochemistry 1979; 18(5):899-904.
- Marlar RA, Kleiss AJ, Griffin JH. Human Protein C: inactivation of factor V and VIII in plasma by the activated molecule. Ann NY Acad Sci. 1981; 370:303-310.
- Comp PC, Esmon CT. Recurrent venous thromboembolism in patients with a partial deficiency of Protein S. New Engl J Med 1984; 311(24):1525-1528.
- Schwarz HP, Fischer M, Hopmeier P, Batard MA, Griffin JH. Plasma protein S deficiency in familial thrombotic disease. Blood 1984; 64(6):1297-1300.
- Kosch A, Junker R, Wermes C, Nowak-Gottl U. Recurrent pulmonary embolism in a 13-year-old male homozygous for the prothrombin G20210A mutation combined with protein S deficiency and increased lipoprotein (a). Thromb Res 2002; 105(1):49-53.
- D’Angelo A, Vigano-D’Angelo S, Esmon CT, Comp PC. Acquired deficiencies of protein S. Protein S activity during oral anticoagulation, in liver disease, and in disseminated intravascular coagulation. J Clin Invest 1988; 81:1445-1454.
- Boerger LM, Morris PC, Thurnau GR, Esmon CT, Comp PC. Oral contraceptives and gender affect Protein S status. Blood 1987; 69(2):692-694.
- Melissari E, Kakkar VV. The effects of oestrogen administration on the plasma free protein S and C4b-binding protein. Thromb Res 1988; 49(5):489-495.
- Garcia de Frutos P, Alim RI, Hardig Y, Zoller B, Dählback B. Differential regulation of alpha and beta chains of C4b-binding protein during acute-phase response resulting in stable plasma levels of free anticoagulant protein S. Blood 1994; 84(3):815-822.
- Comp PC, Thurnau GR, Welsh J, Esmon CT. Functional and immunologic protein S levels are decreased during pregnancy. Blood 1986; 68(4):881-885.
- Fernandez JA, Estelles A, Gilabert J, Espana F, Aznar J. Functional and immunologic protein S in normal pregnant women and in full-term newborns. Thromb Haemost 1989; 61(3):474-478.
- Kisiel W, Hermodson MA, Davie EW. Factor X activating enzyme from Russell’s viper venom: isolation and characterization. Biochemistry 1976; 15(22):4901-4906.
- Biosafety in Microbiological and Biomedical Laboratories 5th ed. Centers for Disease Control and Prevention / National Institutes of Health, 2009.
- Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guidelines – Fifth Edition, CLSI, M21-A5; Vol. 28, No. 5, 2008.
- Cembrowski GS, Carey RN. Laboratory Quality Management. ASCP Press 1989; 166-172.
- CLIA 2004 – Code of Federal Regulations, 42CFR493.1269, 2004.
- Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline-Third Edition, CLSI, C28-A3; Vol. 28, No. 30, 2008.

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